



Your partner in saving sight.®

## FitKit® Mini-Guide

**FOR 16-19mm SCLERAL LENSES**

BostonSight®  
464 Hillside Avenue, Suite 205  
Needham, MA 02494

888-SCLERAL (888-725-3725)  
[BostonSightSCLERAL.org](http://BostonSightSCLERAL.org)

# 18, 18.5 and 19mm FitKit®

19.0 mm  
3.2 Sag  
E-R90-1e1  
FLAT

19.0 mm  
3.2 Sag  
E-R90-2e1  
STANDARD

19.0 mm  
3.2 Sag  
E-R90-3e1  
STEEP

18.5 mm  
3.0 Sag  
E-R85-4e1  
FLAT

18.5 mm  
3.0 Sag  
E-R85-5e1  
STANDARD

18.5 mm  
3.0 Sag  
E-R85-6e1  
STEEP

18.0 mm  
2.8 Sag  
E-R80-7e1  
FLAT

18.0 mm  
2.8 Sag  
E-R80-8e1  
STANDARD

18.0 mm  
2.8 Sag  
E-R80-9e1  
STEEP

19.0 mm  
3.6 Sag  
E-R90-10e1  
STANDARD

SmartSight®  
Lenses  
  
Additional front surface  
eccentricity options for  
refractive purposes

18.5 mm  
3.4 Sag  
E-R85-11e1  
STANDARD

18.0 mm  
3.2 Sag  
E-R80-12e1  
STANDARD

18.0 mm  
3.2 Sag  
FSE0  
E-R80-12e0  
STANDARD

SmartSight®  
Lenses  
  
Additional front surface  
eccentricity options for  
refractive purposes

FITTING LENSES  
RIGHT

All trial lenses have a power of +0.25D, base curve of 8mm, built-in FSE1, and built-in oval optic zones.

Nine core lenses

Additional sag lenses

Additional front surface eccentricity (FSE) lenses for re-fractive purposes.

19.0 mm  
3.6 Sag  
E-L90-10e1  
STANDARD

19.0 mm  
3.2 Sag  
E-L90-1e1  
FLAT

19.0 mm  
3.2 Sag  
E-L90-2e1  
STANDARD

19.0 mm  
3.2 Sag  
E-L90-3e1  
STEEP

18.5 mm  
3.4 Sag  
E-L85-11e1  
STANDARD

18.5 mm  
3.0 Sag  
E-L85-4e1  
FLAT

18.5 mm  
3.0 Sag  
E-L85-5e1  
STANDARD

18.5 mm  
3.0 Sag  
E-L85-6e1  
STEEP

18.0 mm  
3.2 Sag  
FSE2  
E-L80-12e2

18.0 mm  
3.2 Sag  
E-L80-12e1  
STANDARD

18.0 mm  
2.8 Sag  
E-L80-7e1  
FLAT

18.0 mm  
2.8 Sag  
E-L80-8e1  
STANDARD

18.0 mm  
2.8 Sag  
E-L80-9e1  
STEEP

SmartSight®  
Lenses  
  
Additional front surface  
eccentricity options for  
refractive purposes

FITTING LENSES  
LEFT

2 | FitKit Mini-Guide

# 16, 16.5, and 17mm FitKit®

## FITTING LENSES RIGHT

All trial lenses have a power of +0.25D, base curve of 8mm, built-in FSE1, and built-in oval optic zones.

<b>17.0 mm</b> 3.0 Sag D-R70- <b>1e1</b> FLAT	<b>17.0 mm</b> 3.0 Sag D-R70- <b>2e1</b> STANDARD	<b>17.0 mm</b> 3.0 Sag D-R70- <b>3e1</b> STEEP	<b>17.0 mm</b> 3.4 Sag D-R70- <b>10e1</b> STANDARD	<b>SmartSight® Lenses</b>  Additional front surface eccentricity options for refractive purposes
<b>16.5 mm</b> 2.8 Sag D-R65- <b>4e1</b> FLAT	<b>16.5 mm</b> 2.8 Sag D-R65- <b>5e1</b> STANDARD	<b>16.5 mm</b> 2.8 Sag D-R65- <b>6e1</b> STEEP	<b>16.5 mm</b> 3.4 Sag D-R65- <b>11e1</b> STANDARD	<b>16.0 mm</b> 3.0 Sag FSE0 D-R60- <b>12e0</b> STANDARD
<b>16.0 mm</b> 2.6 Sag D-R60- <b>7e1</b> FLAT	<b>16.0 mm</b> 2.6 Sag D-R60- <b>8e1</b> STANDARD	<b>16.0 mm</b> 2.6 Sag D-R60- <b>9e1</b> STEEP	<b>16.0 mm</b> 3.0 Sag D-R60- <b>12e1</b> STANDARD	<b>16.0 mm</b> 3.0 Sag FSE2 D-R60- <b>12e2</b> STANDARD

## FITTING LENSES LEFT

<b>SmartSight® Lenses</b>  Additional front surface eccentricity options for refractive purposes	<b>17.0 mm</b> 3.4 Sag D-L70- <b>10e1</b> STANDARD	<b>17.0 mm</b> 3.0 Sag D-L70- <b>1e1</b> FLAT	<b>17.0 mm</b> 3.0 Sag D-L70- <b>2e1</b> STANDARD	<b>17.0 mm</b> 3.0 Sag D-L70- <b>3e1</b> STEEP
<b>16.0 mm</b> 3.0 Sag FSE0 D-L60- <b>12e0</b>	<b>16.5 mm</b> 3.2 Sag D-L65- <b>11e1</b> STANDARD	<b>16.5 mm</b> 2.8 Sag D-L65- <b>4e1</b> FLAT	<b>16.5 mm</b> 2.8 Sag D-L65- <b>5e1</b> STANDARD	<b>16.5 mm</b> 2.8 Sag D-L65- <b>6e1</b> STEEP
<b>16.0 mm</b> 3.0 Sag FSE2 D-L60- <b>12e2</b>	<b>16.0 mm</b> 3.0 Sag D-L60- <b>7e1</b> STANDARD	<b>16.0 mm</b> 2.6 Sag D-L60- <b>7e1</b> FLAT	<b>16.0 mm</b> 2.6 Sag D-L60- <b>8e1</b> STANDARD	<b>16.0 mm</b> 2.6 Sag D-L60- <b>9e1</b> STEEP

# Lens Handling and Conditioning

Trial lenses in the FitKit will be shipped wet to ensure proper wetting. The conditioning solution must be replaced with fresh solution (i.e. Boston Simplus) upon FitKit receipt and must be replaced every 30 days. Proper cleaning, disinfection and handling is required before each use.


- ① Wash hands thoroughly with non-moisturizing soap and water before handling lenses.
- ② Place one lens in the palm of your hand. Place one or two drops of your daily cleaning solution on the outside surface of the lens and rub gently for 15 seconds. Flip the lens over and do the same on the inside surface of the lens.
- ③ Rinse the lens thoroughly with preservative-free saline.
- ④ Place lens in conditioning solution, such as Boston Simplus, for at least 1 hour before use.

Your BostonSight SCLERAL trial lens is now ready for patient application.

## Fitting Algorithm Schematic

- ① Choose diameter.
- ② Start at the center with a standard shape lens and standard sagittal depth at each given diameter.
- ③ If the standard shape lens fits loose, choose the steep shape lens.
- ④ If the standard shape lens fits tight, choose the flat shape lens.

### Fitting guide based on a simplified system for each eye



## BETTER WEIGHT DISTRIBUTION FOR HIGH SAG CORNEAS

### 18mm

- Normal cornea with ocular surface disease
- Post surgical cornea with moderate ectasia
- Compromised/fragile grafts
- Chronic exposure
- Moderate keratoconus
- Limbal stem cell deficiency
- HVID > 11.50mm

### 18.5 – 19mm

- Advanced keratoconus (Approx. 65D+)
- Post surgical cornea with highly irregular ectasia
- Advanced ocular surface disease
  - GVHD
  - SJS
  - Sjögren's

# Fitting Principles

The process of fitting the BostonSight SCLERAL lens is based on identifying the best fitting trial lenses and adapting their geometries and power to create eye-specific lenses.

- ① Identify the trial lens having the best initial fit using the fitting algorithm.
- ② Re-evaluate the fit after the appropriate settling time.
- ③ When indicated, replace with one having more appropriate parameters.
- ④ Repeat the process until the best fitting trial lens is identified.
- ⑤ Perform spherical over-refraction to determine lens power. Vision may be optimized using a different trial with different front eccentricity values or by performing sphero-cylindrical refraction over the final diagnostic lens.
- ⑥ Order the lens with any modifications as needed in FitConnect®.

# Lens Identification

Easy diagnostic lens ID with unique laser-engraved model number

**E-R 80-1 e1**

LENS SET  
SERIES

EYE  
R VS. L

DIAMETER

LENS NUM-  
BER 1-2

ECCENTRICITY  
FSE0, FSE1 OR  
FSE2

Example Shown:

Lens Set Series: E

Right Lens

18.0mm diameter

Lens #1

FSE1



Laser engraved model  
number location



Your partner in saving sight.\*

BostonSight®  
464 Hillside Avenue, Suite 205  
Needham, MA 02494

Copyright © 2025

888-SCLERAL (888-725-3725)  
BostonSightSCLERAL.org

MKT-020 Rev 1