



BOSTONSIGHT

SCLERAL

Going Beyond the Fit™

with our New Data-Driven
Lens Fitting System



INNOVATION | EDUCATION

Introducing BostonSight® SCLERAL

Data-Driven for an Improved Fit

6 YEARS 

~7,000 EYES



Welcome to a new vision for irregular corneas and ocular surface disease

BostonSight SCLERAL is the **first** scleral lens fitting system based on scleral anatomy, driven by clinical data, and the first with front surface eccentricity options to improve both quantitative and qualitative vision.

Our Goal: Provide your patients comfort and improved ocular health with an accurate and dependable fit through:

- Minimized modifications
- Reduced number of patient sessions

Addressing a Spectrum of Need

BostonSight SCLERAL is appropriate for a wide range of conditions including:



IRREGULAR OR DISTORTED CORNEA



OCULAR SURFACE DISEASE



AFTER SURGERY

A More Dependable Fit, Right from the Start

Each fitting uses a single starting point. BostonSight SCLERAL is the first lens fitting system with right and left eye anatomical designs. The fitter customizes with options including:



Quadrant-specific toric peripheral haptic system (PHS™)



Front surface eccentricity to improve both quantitative and qualitative vision



18.0mm, 18.5mm, and 19.0mm diameter designs

Fewer fitting modifications help ensure a predictable and accurate outcome for your patients. What's more, less chair time per patient also helps make the overall fitting process more efficient.



Going Beyond the Fit™

Comprehensive Support for an Enhanced Experience

BostonSight SCLERAL's Beyond the Fit program includes:



FitConnect™: An intuitive, web-based fitting and order management system for BostonSight SCLERAL which empowers the practitioner to:

- Make modifications in real-time for more accurate endpoints
- Adjust designs using graphical representations of individual lens profiles
- Modify sagittal height and 4 independent hemi-meridians



BostonSight SCLERAL Vouchers: BostonSight is reducing financial barriers to BostonSight SCLERAL by offering a lens voucher program



MyFitReport™: Provides patient-generated health data (PGHD), with patient-specific feedback, to allow practitioners to evaluate product performance and patient satisfaction

Focus on what Matters: Treating Patients

Delivering a system with a detailed level of customization that optimizes fitting efficiency and increases satisfaction for both you and your patients:

- ✓ No start-up costs required
- ✓ No extensive training needed
- ✓ No special equipment necessary
- ✓ Educational materials provided for ocular health and scleral lenses

BostonSight SCLERAL Lens Parameter Availability

Diameters	18.0mm, 18.5mm, and 19.0mm
Sphere power	-20.00 Diopters to +20.00 Diopters
Sagittal height	2.0mm to 6.0mm in 0.1mm (100 μ m) steps
Peripheral Haptic System (PHS™)	Customizable
Polymer material	Boston® EQII Boston® XO2 Optimum-Extra Optimum-Extreme
Options	Quadrant-specific toric PHS™ <ul style="list-style-type: none">• Front surface eccentricity• Front surface toric Rx



For more information on the BostonSight SCLERAL fitting system, to set-up an account, or to request a live FitConnect demo, contact us at:

www.BostonSightSCLERAL.org

888-SCLERAL



BOSTONSIGHT

BostonSight, a nonprofit organization, is known worldwide as a pioneer in the field of scleral lenses and as the innovator behind life-changing BostonSight PROSE treatment. BostonSight is dedicated to improving the quality of life for patients and their families by transforming the understanding, treatment, and care of the ocular surface.